

# Use of the DEPArray platform to detect, isolate, and molecularly characterize pure tumor cells from peripheral blood samples enriched using the CellSearch® system

G Medoro<sup>1</sup>, S Gross<sup>2</sup>, N Manaresi<sup>1</sup>, M Sergio<sup>1</sup>, F Fontana<sup>1</sup>, S Gianni<sup>1</sup>, A Calanca<sup>1</sup>, E Peruzzi<sup>1</sup>, M Banzi<sup>1</sup>, G Signorini<sup>1</sup>, C Rao<sup>2</sup>, J Patel<sup>3</sup>, J Karkera<sup>3</sup>, G Giorgini<sup>1</sup>, M Mata<sup>3</sup> and M C Connelly<sup>2</sup>

Abstract #10616

<sup>1</sup>Silicon Biosystems Spa, Bologna, Italy <sup>2</sup>Veridex, LLC, Huntingdon Valley, PA <sup>3</sup>Centocor, Radnor, PA

## Abstract

**Background:** Circulating tumor cells (CTC) offer the potential for serially monitoring the molecular profile of a tumor. However, enrichment techniques provide a level of purity problematic for most molecular analysis methods, and do not readily provide for analysis of tumor cell heterogeneity. We evaluate the use of DEPArray™ (Silicon Biosystems), an automated system enabling image-based cell sorting with single-cell resolution, for CTC isolation and characterization from enriched blood samples.

**Methods:** Experiments were carried out with healthy-donor blood (HB) collected in CellSave tubes, spiked with tumor cells (TC) and enriched on Veridex's AutoPrep with the CellSearch® Epithelial Cell Kit. DEPArray™ system was used for detection and multiple recoveries of single TCs (or control WBCs) and 5 cell batches. A comparison of blind enumeration results with Veridex's CellTracks Analyzer II® (CTAIL) was carried out on replicate samples.

**Results:** No TCs were detected among negative controls (n=10). TC count (normalized to sample volume analyzed) was compared sample-wise for each replicate (n=20): DEPArray/CTAIL count was, on average, 100% (standard deviation = 52%). Enriched mixtures of Her2+ and Her2- TCs spiked in HB samples (n=5), were sorted by DEPArray and recovered into separate tubes. By phenotypical re-analysis no Her2+ cells were detected among the Her2- cell fraction and vice versa, neither were donor WBCs found (100% purity). KRAS-mutated, A549 cells spiked in HB were enriched and loaded on DEPArray. Individual fractions containing either 1 to 5 tumor cells or donor WBCs were sorted. Whole Genome Amplification (Ampli1™ WGA, Silicon Biosystems), KRAS specific gene amplification and Capillary Electrophoresis sequencing were carried out. TCs successfully amplified showed only mutated KRAS, (WBCs were only wild-type).

**Conclusions:** DEPArray™ achieved 100% purity, eliminating all white blood cells (WBC), in the isolation of a mixed population of tumor cell lines downstream of CellSearch® enrichment. This enabled molecular profiling of pure tumor cells from whole blood spiked tumor cell lines. Detection of molecular heterogeneity of tumor cells is demonstrated through KRAS sequencing

## Tumor Cell Detection

**Experiment work-flow**

**DEPArray™ Cell Browser screenshot**

**Image gallery of TCs from DEPArray™**

**Summary of results**

Cells prepared in Veridex central lab, PA (V-HV); count performed in blind by remote labs: Veridex Enschede, NL (V-EU) and Silicon Biosystems (SB) with DEPArray. Negative controls (only Healthy donor Blood, n=10), were negative in all platforms (data not shown). For spiking experiments (n=20, Table 1) of TCs (SKBr3), a normalization coefficient (9,26/15=61,7%) is used on Veridex counts to take into account volume effectively analyzed by DEPArray (Main Chamber volume). DEPArray/CTAIL remote count Ratio = SB/(V-EU)=100%, σ=52%. Average number of TC in remote labs [(SB+V-EU)/2]/(V-HV)=80% due to shipping loss, σ=16%, (<σ(SB/V-EU) due to negative correlation from sample splitting after pooling).

**Cells (Normalized)**

V-HV	V-EU	SB	SB/EU	(EU+SB)/2
12	12	12	97%	99%
12	9	13	150%	88%
94	91	25	27%	62%
10	3	6	194%	43%
12	16	9	56%	107%
9	7	5	74%	64%
117	115	61	53%	75%
132	105	78	74%	69%
93	83	45	54%	69%
69	65	60	93%	90%
102	103	85	82%	92%
48	27	55	202%	85%
159	146	142	97%	90%
162	71	150	211%	68%
66	59	48	82%	81%
44	41	38	92%	91%
44	31	22	71%	59%
62	70	46	66%	93%
78	57	57	100%	73%
80	73	90	123%	102%

**Table 1**

mean	100%	80%
std-dev	52%	16%

## Cell Separation

**Experiment work-flow**

**Separation by Immunophenotype: Results evaluation by CTAIL®**

Duplicate preparations of Her2+/Her2- admixed cells (SKBr3/PC3-9) spiked in HB (n=5) were enriched with CellSearch® and scanned by CTAIL® before and after DEPArray sorting of Her2+ and Her2- fractions. Sorted fractions re-injected in Veridex cartridges and analyzed with CTAIL revealed 100% purity, with no WBCs or wrong type of tumor cells

**CellTracks Analyzer II® Images**

**Analyzed by CTAIL (V-HV)**

Replicate	Sorted by DEPArray™ (SB)	Her2-	Her2+	Purity	
1	Her2-	50	28	0	100%
	Her2+	47	0	17	100%
2	Her2-	6	3	0	100%
	Her2+	27	0	16	100%
3	Her2-	17	7	0	100%
	Her2+	26	0	16	100%
4	Her2-	5	1	0	100%
	Her2+	11	0	6	100%
5	Her2-	8	5	0	100%
	Her2+	16	0	11	100%

## Cell Molecular Characterization

**Experiment work-flow**

**Molecular characterization for purity and mutational status**

Across four Spiking experiments of viable (>90%) KRAS mutated cell lines (SW480 n=2, A549 n=1, mix SW480/A549 n=1) in HB, multiple recoveries (range 12-21 per experiment) of individual cells (n=56), 5 cell batches (n=5), or negative controls (n=8) were carried out. Ampli1™ WGA kit (Silicon Biosystems) products from each tube were DNA-fingerprinted (home-brewn 11 loci multiplex reaction, to confirm cell presence, identity and purity) and KRAS gene-specific amplification products were sequenced (on ABI 3730xl). In 91% (51/56) of single cell recoveries, cell presence was confirmed. The 5 tubes with no signal from STR and KRAS suggests that the cell has been removed during supernatant removal before WGA. All successfully amplified cells matched 100% KRAS mutational status and DNA fingerprint, no signals was detected in negative controls recoveries (buffer only). In the mixed tumor cells experiment, different KRAS mutations (and DNA fingerprints) were detected in different tumor cell recoveries, reflecting cell heterogeneity.

**Molecular Analysis Results**

**Summary**

	N	KRAS call	KRAS no call	KRAS call rate	# Alleles expected	# Alleles found	ACR % (non void)
Single cells	56	51	5	91%	957	588	68.6%
TC	38	36	2	95%	557	388	69.7%
WBC	18	15	3	83%	300	200	66.7%
5 cells	5	5	0	100%	97	82	84.6%
TC	1	1	0	100%	17	15	88.2%
WBC	4	4	0	100%	80	67	83.8%
Negative Controls	8	0	0	NA	0	0	NA

**KRAS sequencing**

**DNA fingerprinting**

**Sample Classification Table**

Sample ID	Recovery ID	# cells	cell type	KRAS PCR	KRAS type	Notes	# Alleles expected	# Alleles found	Allele Call Rate %	Result Classification
VRX206R2-SW480	1	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	2	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	3	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	4	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	5	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	6	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	7	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	8	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	9	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	10	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	11	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	12	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
VRX206R3-A549	1	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	2	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	3	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	4	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	5	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	6	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	7	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	8	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	9	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	10	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	11	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
	12	1	TC	1	G12V	allelic count refer to SW480	12	13	100%	MATCH
VRX206R3-WBC	1	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	2	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	3	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	4	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	5	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	6	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	7	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	8	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	9	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	10	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	11	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL
	12	1	WBC	1	WT	allelic count refer to WBC	20	14	70%	NO-CALL